

AUTOMATIC PRESSURE INFUSOR AUTOPRESS

PRECISE PRESSURE SETTING. PRESET PRESSURE.

BIEGLER **Autopress** is another component in our range of warming equipment. It can be used either directly with our blood and infusion warmers or as a stand-alone solution wherever fluids have to be dispensed under constant pressure.



AUTOMATIC PRESSURE INFUSOR AUTOPRESS

- Preset pressure is always maintained
- Precise pressure setting due to electronic control
- Accomodated pressure cuffs 2 x 500 or 2 x 1000 cc
- Fits on any IV pole
- Small and lightweight
- Pressure limited to 300 mmHg
- Mains or battery operated
- Pressure cuffs can be emptied rapidly and easily
- Significant reduction of set up time
- Battery backup up to 90 minutes

TECHNICAL DATA:

Device:	automatic pressure infusor
Type designation:	Autopress
Voltage:	9V DC
Power consumption:	max. 12 W
Type:	mains and rechargeable battery operated
Classification:	Ila according to rule 11
Degree of protection against electric shock:	B
Dimension:	90 x 115 x 260 mm
Weight:	1.46 kg
Type of operation:	continuous operation
Mains adaptor Type designation:	12W P2EFMW3
Input:	110-240V AC, 50/60 Hz, 250 mA
Output:	9V DC, 1.33A max.
Accessories:	Pressure cuff 500 cc Pressure cuff 1000 cc

AUTOPRESS

BIEGLER GmbH
Allhangstrasse 18a
3001 Mauerbach, Austria

Tel: +43 1 979 21 05
Fax: +43 1 979 21 05 16
office@biegler.com
www.biegler.com

BIEGLER
MEDIZIN ELEKTRONIK

CE
0123

