

EMERGENCY WARMER ESH 04

VERSATILE. FOR MOBILE USE.

The Biegler **ESH 04** emergency warmer has been specifically designed for mobile use. It can keep two pre-warmed 500 ccm infusion bottles at a control temperature of 37°C at the place of use for up to 90 minutes. The device is suitable for continuous operation when mains-connected, and it can also pre-warm infusion bottles in this mode.



EMERGENCY WARMER ESH 04

- Designed for mobile use
- Versatile, compact size, low weight
- Fits into emergency kit
- Maintains the temperature for up to 90 minutes at the site of use
- 24 h continuous operation possible with mains power
- Integrated pressure cuff
- Solution dispensed directly from warmer

TECHNICAL DATA:

| | |
|---|---|
| Device: | Emergency warmer |
| Type designation: | ESH 04 |
| Voltage: | 12V DC |
| Power consumption: | max. 20 W |
| Type: | mains and rechargeable battery operated |
| Degree of protection against ingress of liquids: | IPX splash proof |
| Degree of protection against electric shock: | B |
| Dimensions: | 160 x 105 x 300 mm |
| Weight: | 1.3 kg |
| Classification: | IIa according to rule 9 |
| Type of operation: | continuous operation |

| | |
|--------------------------|----------------------------------|
| Mains adaptor | |
| Type designation: | 24W P2EFSW3 |
| Input: | 110-240V AC, 50-60 Hz, 550 mA |
| Output: | 12V DC, 2A max. |

ESH 04

BIEGLER GmbH
Allhangstrasse 18a
3001 Mauerbach, Austria

Tel: +43 1 979 21 05
Fax: +43 1 979 21 05 16
office@biegler.com
www.biegler.com

BIEGLER
MEDIZIN ELEKTRONIK

CE
0123

