

# **AUTOMATIC PRESSURE INFUSOR**

## **PRECISE PRESSURE SETTING.**

## **PRESET PRESSURE.**

*BIEGLER AUTOPRESS* is another component in our range of warming equipment. It can be used either directly with our blood and infusion warmers or as a stand-alone solution wherever fluids have to be dispensed under constant pressure.



## AUTOMATIC PRESSURE INFUSOR



- Preset pressure is always maintained
- Precise pressure setting due to electronic control
- Accomodates pressure cuffs 2 x 500, 2 x 1000 or 2 x 3000 ml
- Fits on any IV pole and rail
- Small and lightweight
- Pressure limited to 300 mmHg
- autopress: Mains operated  
autopress >e: Mains operated with up to 3h battery backup
- Pressure cuffs can be emptied rapidly and easily
- Significant reduction of set up time

### TECHNICAL DATA

<b>Device:</b>	automatic pressure infusor	
<b>Type designation:</b>	autopress	autopress >e
<b>Power consumption:</b>	max. 36 VA	
<b>Operating voltage:</b>	110 V-240 V AC, 50-60 Hz	
<b>Supply type:</b>	mains operated	mains and battery operated
<b>Classification (93/42/EEC):</b>	IIa according to Rule 9	
<b>Degree of protection against electric shock:</b>	Type CF defibrillationproof	
<b>Operation mode:</b>	continuous operation	
<b>Dimension:</b>	100 x 230 x 180 mm	
<b>Weight:</b>	2.1 kg	2.3 kg
<b>Accessories:</b>	Pressure cuff 500 ml Pressure cuff 1000 ml Pressure cuff 3000 ml	