

EMERGENCY WARMER ESH 04

VERSATILE. FOR MOBILE USE.

The Biegler **ESH 04** emergency warmer has been specifically designed for mobile use. It can keep two pre-warmed 500 ccm infusion bottles at a control temperature of 37°C at the place of use for up to 90 minutes. The device is suitable for continuous operation when mains-connected, and it can also pre-warm infusion bottles in this mode.





ESH 04



EMERGENCY WARMER ESH 04

- Designed for mobile use
- Versatile, compact size, low weight
- Fits into emergency kit
- Maintains the temperature for up to 90 minutes at the site of use
- 24 h continuous operation possible with mains power
- Integrated pressure cuff
- Solution dispensed directly from warmer

TECHNICAL DATA:

Device:	Emergency warmer
Type designation:	ESH 04
Voltage:	12V DC
Power consumption:	max. 20 W
Type:	mains and rechargeable battery operated
Degree of protection against ingress of liquids:	IPX splash proof
Degree of protection against electric shock:	B
Dimensions:	160 x 105 x 300 mm
Weight:	1.3 kg
Classification:	Ila according to rule 9
Type of operation:	continuous operation
Mains adaptor	
Type designation:	24W P2EFSW3
Input:	110-240V AC, 50-60 Hz, 550 mA
Output:	12V DC, 2A max.

BIEGLER GmbH
Allhangstrasse 18a
3001 Mauerbach, Austria

Tel: +43 1 979 21 05
Fax: +43 1 979 21 05 16
office@biegler.com
www.biegler.com

BIEGLER
MEDIZIN ELEKTRONIK

CE
0123

